



Standard Operating Procedure (SOP)

Filmerics F40

(MET-04)

In case of emergency please call 911

For any other major safety concern contact EHRS at: 215-898-4453 or via email: ehrs@ehrs.upenn.edu

If there is an error on the system/tool please report it on IRIS, we will take care of it

Please *DO NOT* run diagnosis without a staff member's approval

General safety tips and common mistakes

1) Make sure that you are logged into the tool on IRIS before use. If you do not log in, you cannot run the tool.

Filmetrics F40



Procedure Overview

- 1 Check-in
- 2 Start-up

3 Baseline and reference check

3.1 Baseline setting

3.1.1 Measure Tab

3.1.2 Focus the camera/optics

3.2 Measurement on the standard SiO2 sample

- 3.2.1 Step 1
- 3.2.2 Step 2
- 3.2.3 Step 3
- 3.2.4 Measurement

4 Measurement on Sample

- 4.1 Step 1
- 4.2 Step 2
- 4.3 Step 3
- 4.4 Creating/editing measurement recipes
- 4.5 Measuring your sample and reviewing data
- 5 Leaving the tool in an "idle" state

6 Check-out



3. Baseline and reference
check
To obtain accurate
measurements, baseline
and reference check the
F40 at the beginning of
each tool reservation
and/or when switching
between thin film types
(e.g.: Oxide, then
photoresist)







3.2.2 Step 2

3.2.2.1. The dialog box for step 2 is opened when step 1 is completed.





3.2.2.3. Focus on the
"Reference" field using the
CCD image in the "Live
Video" tab, if necessary.
3.2.2.4. Choose "Si" from
the pull-down menu of the
"Reflectance Standard" box.
3.2.2.5. Click the "Take
Reflectance Standard"
button.

BK7
Si
Enter Constant
Acrylic
AI2O3
Au
BK7
BSG
CaF2
Fused Silica (2-Sided)
GaAs
GaAs_hot
Ge
InP
NiP
Polycarbonate
Polyolefin
Polystyrene
Quartz
Si
SiO2
Clear Most Recent List







4.1.3. Place your sample beneath the objective lens.
4.1.4. Focus on the sample using the CCD image in the "Live Video" tab, if necessary.
4.1.5. Click the "Take Sample Reflectance" button.







4.4 Creating/editing ame: SiO2 on Si • Author: Last Modified: 9/25/2013 4:42:15 PM ilm Stack Analysis Options Alarms Ac measurement recipes Units: Angstroms (Å) 💌 Medium Air . 4.4.1. Select desired recipe 1 SiO2 rate Si from pull-down menu above the "Edit Recipe" button and click "OK". 4.4.2. Click "Edit Recipe" to open the "Edit Recipe" OK Cancel Apply dialog box. 4.4.3. Click "film stack" to Note: To change or add the parameters, check Meas box, and change the add, remove, or change film parameters. types in the stack to be measured Edit Recipe 4.4.4. When done with 💕 🎦 🔒 🔣 🕅 edits, click "apply" then Recipe Name: SiO2 on Si "OK" to apply changes Save Recipe as different filename temporarily Film Stack Analysis Options Alarms Wafer Ma 4.4.5. To save an edited recipe as a new recipe, click "apply", then "save recipe as different filename" (represented by the "disk with pencil" icon), rename, then click "OK"

4.5 Measuring your sample and reviewing data

4.5.1. Place your sample beneath the objective lens.4.5.2. Focus on your sample.

4.5.3. Select the desired measurement recipe, edit and/or "save as", if desired 4.5.4. Use live video feed and X & Y micrometers to locate area to be measured

Use "center" to move spot measurement location to center of live video display

Use "settings" to adjust spot size for measuring, as needed

Use "zoom" to help locate smaller features/areas to be measured

4.5.5. Click "Measure" and view results in the lower right hand corner of the screen

"idle" state 5.1. When done using the F40, turn off the light

5. Leaving the tool in an

source only, unless another





will be using the tool immediately after you 5.2. Leave the FILMeasure software open 5.3. Leave main unit power on	
6. Check-out	
Log-off the IRIS scheduler	

Feel free to contact the staff members with any questions about your process and the tool.

Last modified: 12/6/2021 by Hiromichi Yamamoto