



11-9-2015

## PDMS-PDMS Bonding Protocol - Technics

Steven Henry  
*University of Pennsylvania*

Follow this and additional works at: [http://repository.upenn.edu/scn\\_protocols](http://repository.upenn.edu/scn_protocols)

---

Henry, Steven, "PDMS-PDMS Bonding Protocol - Technics", *Protocols and Reports*. Paper 3.  
[http://repository.upenn.edu/scn\\_protocols/3](http://repository.upenn.edu/scn_protocols/3)

This paper is posted at ScholarlyCommons. [http://repository.upenn.edu/scn\\_protocols/3](http://repository.upenn.edu/scn_protocols/3)  
For more information, please contact [libraryrepository@pobox.upenn.edu](mailto:libraryrepository@pobox.upenn.edu).

---

# PDMS-PDMS Bonding Protocol - Technics

## **Keywords**

PDMS, Bonding, Technics

## **Creative Commons License**



This work is licensed under a [Creative Commons Attribution 4.0 License](https://creativecommons.org/licenses/by/4.0/).

## PDMS-PDMS Bonding Protocol – Technics

Updated on 11/09/2015

**Critical factors:**

- Set O<sub>2</sub> pressure to 2.16 Torr
- Power to 30W
- Plasma for 30 s
- Wait 20 minutes before testing bond

**Table of Contents:**


1. Materials
2. Process Parameters
3. Protocol
4. Tested Parameter Results

**Materials:**

- 10:1 PDMS spun glass microscope slides:
  - 2.5 min plasma clean per side (2.16 Torr, 30 W, O<sub>2</sub> flow bob at “97.5”)
  - 10g base + 1g cure mix vigorously and degas 1 hour
  - Spin glass coverslides, 10 s @ 500 rpm, 50 s @ 3000 rpm
  - Allow PDMS to level at RT for 2 hrs
  - Cure overnight at 65C
- 20:1 PDMS disks (~10 mm tall, ~12 mm diameter)
  - 10 min sonication in 200 proof EtOH
  - 2X 60 mL rinse in MilliQ water
  - N<sub>2</sub>(g) dried

**Process Parameters:**

- O<sub>2</sub> flow rate value “97.5” (bob between “95” and “100” markings)
- Vacuum (process valve) adjusted until pressure as reported on black box to left of Technics reads 2.16 Torr ± 0.02 Torr
- Power set to 30 W (0.030 kW)
- Plasma distributing plate in place
- Positioning of samples in chamber as follows:



**Protocol:**

- Prior to loading samples, perform run on an empty vacuum chamber to properly set power supply knob position for desired process parameters.
- Load samples in orientation shown above.
- Reestablish process parameters by *manipulation of process valve only*.
- Perform plasma exposure for specified time.
- Vent chamber.
- Place PDMS in conformal contact with PDMS spun glass.
- Apply gentle but uniform pressure to PDMS for 10s.
  - Note: Extremely thin glass coverslips (#1-#2) will warp if excessive pressure is applied. Gentle and uniformly distributed pressure is only necessary to ensure conformal contact of the PDMS/PDMS interface.
- Incubate PDMS on PDMS spun glass at room temperature for 15 min.
  - Note: If more than one device was bonded, do not leave them in contact with one another at least for a few hours as partial plasma activation could have occurred on the other surfaces of the PDMS.

**Tested Parameter Results:**

Power (W)	Pressure (Torr)	O <sub>2</sub> Flow Rate "Value"	Treatment Duration (s)	Incubation Duration (min)	Peel Test Outcome
0	2.16	97.5	10 (held at 2.16 Torr w/o plasma)	15	Fail
30	2.16	97.5	5	15	Fail
30	2.16	97.5	10	15	Fail
30	2.16	97.5	30	15	Pass
30	2.16	97.5	120	15	Pass
30	2.16	97.5	300	15	Pass
30	0.70	47.5	5	15	Pass
30	0.70	47.5	10	15	Pass
30	0.70	47.5	120	15	Fail
60	2.16	97.5	10	15	Pass